

- Variable sizes of staple reload to offer you more options for different tissue thickness.
- Build-in knife blade for each staple reload to provide the complete tissue cutting and reduce the possibility of implantation metastasis after tumour-reductive surgery.
- 6 rows of staple lines to improve the safety of anastomosis.
- Three terraced staple heights into each purple and black cartridge to provide better tissue squeezing force but not over squeeze tissue. The better blood circulation will help anastomosis recover faster.



TCS Series Disposable Endoscopic Cutting Stapler

NO.	Description	Model	Remark
1	TCS series Disposable endoscopic cutting stapler	TCS60	(Short)
2		TCS160	(Regular)
3		TCS260	(Long)

Applicable Staple Reload

TRR Series Disposable Endoscopic Cutting Stapler Reload (Articulated):

Model	Cartridge Color	Anastomosis length (mm)	Staple height (mm)	Staple quantity (pc)
TRR30W	(White)	31	2.5	48
TRR30B	(Blue)	31	3.5	48
TRR30P	(Purple)	31	3.0-3.5-4.0	48
TRR45W	(White)	43	2.5	66
TRR45B	(Blue)	43	3.5	66
TRR45P	(Purple)	43	3.0-3.5-4.0	66
TRR45H	(Black)	43	4.0-4.5-5.0	66
TRR60B	(Blue)	59	3.5	90
TRR60P	(Purple)	59	3.0-3.5-4.0	90
TRR60H	(Black)	59	4.0-4.5-5.0	90